27.04.2024 1 von 1





Resuscitation Dog Casper, weight: 6.7 kg

Order number: AU-02558

With Casper the resuscitation dog, it is possible to practice mouth to muzzle respiration and heart massage in an animal. The tongue can be pulled out to check for a blockage and the heart rate can be checked in the groin.

With a fire or accidents, apart from people there are often also animals involved. Many people have an emotional bond with their dog or cat. As first-aid assistant, it is important to know the basic resuscitation techniques. These procedures can be practiced on Casper the dog. Casper's lungs can easily be replaced, therefore every student can get a new one. The lungs fall over the muzzle, therefore you do not need to disinfect. A pump is also provided to simulate the heart rate manually.

size: 67.3 x 58 x 23 cm

E-Mail info@rescue-tec.de

Fax.:

+49-6482-608900

+49-6482-608920

weight: 6.7 kg