



## **PAX Emergency Backpack Berlin - Magnet, PAX-Tec (blue, red), PAX-Plan (yellow)**

**Order number:** RD-07515

### **ENOUGH SPACE AND STILL COMPACT**

Our emergency backpack Berlin forms the link between our Wasserkuppe backpacks and the Feldberg range. It is compact and still offers you enough space for good emergency equipment for your daily use. The front pockets are set inside. In them and in the optionally available, colour-coded Velcro pockets, you can store your equipment thematically separated, in accordance with the ABC scheme.

### **FRONT POCKETS, WHICH ARE ALSO INNER POCKETS**

To keep the backpack narrow, we have placed the front pockets inside. You open it from the outside with a double zipper. From the inside you also access the front pocket. If necessary, you can reach into the body of the backpack via the front pocket without opening it completely. Infusions and accessories are stored in the double loop rows of the lower front pocket. Diagnostic materials, for example, fit into the upper front pocket.

### **PLENTY OF RUBBER LOOPS**

In the body you will find numerous rubber loops (which you can take out or exchange for shorter ones). You can attach smaller parts such as respiratory masks quickly and easily. You can stow your gloves, for example, for DeepL access in the tried and tested transparent pockets that are firmly attached to the frame

### **Excellence - EVEN FROM THE OUTSIDE**

Of course we have taken just as much trouble with the exterior design. Your emergency backpack Berlin has a handle on the top and on the side. You can replace the shoulder straps (with quick release buckle). To protect your backpack from moisture, the bottom is made of PAX-Lon. Even if you let j. w. d. he stands on his four firmly riveted feet at a safe distance from dirt and moisture. You can replace the feet yourself if necessary. Or we can do it for you in our repair service. To protect the sliders of the zipper from wear and tear, park them in a zipper garage with the rucksack open.



When closed, this garage protects the zipper chain from abrasion.

The fixation by magnets not only facilitates the removal of the modules, it also promotes easier disinfection.

material: PAX-Tec (blue, red), PAX-Plan (yellow)

dimensions: 57 x 39 x 26 cm

weight: 3.4 kg